

Digital Spectrophotometer



Product Categories: [Analytical Lab Instruments](#)

Product Tags: [Analytical Laboratory Instrument](#)

Product Page: <https://www.labappara.com/product/digital-spectrophotometer-2/>

Product Description

Digital Spectrophotometer - Analytical Laboratory Instrument] - Range
340 960 nm with facility .

The automatic concentration, % .

Transmission and absorbency measurements 3 and half
digit data display.

- separate 3 digit LED display for wave.

length 1 nm resolution with a set of 1 Pair test tubes.

10mm & 50mm path length cuvettes at Extra cost.