

## DIGITAL SPECTROPHOTOMETER



**Product Categories:** [Analytical Lab Instruments](#), [Spectrophotometer](#)

**Product Tags:** [Analytical lab equipment](#), [Spectrophotometer](#)

**Product Page:** <https://www.labappara.com/product/digital-spectrophotometer-3/>

### Product Description Digital Spectrophotometer

Wavelength Resolution being 1 nm.

Wavelength - 340 to 960nm

10 & 50MM path length.

Direct concentration Facility.

Dual Digital display available for wavelength & data.

High durability and strength.

In available wide range.

Standard rates exist.

Comes with optional features of Corvettes of 10mm and 50mm in glass and quartz.

Light Source: Tungsten-Halogen Lamp

Detector: Silicon Photodiode

Resolution: 0.1%T Abs 0.001

Wavelength resolution: 1nm

Accuracy:  $\pm 2.5\text{nm}$

Spectral bandwidth: 5nm