

## Digital Tds Meter



**Product Categories:** [Analytical Lab Instruments](#)

**Product Tags:** [Analytical Laboratory Equipment](#)

**Product Page:** <https://www.labappara.com/product/digital-tds-meter-2/>

### Product Description

Digital Tds Meter - Analytical Laboratory Equipment

Three digit LED display with facility

Digital TDS Meter temperature compensation and digital cell constant adjustment along with TDS cell & other accessories.

### Description

TDS meter

TDS meter

## Measuring range

0~1999ppm

0~19.99ppt

## Accuracy

±2%FS

±2%FS

## Resolution

1ppm

0.01ppt

## Display

4-digital LCD with Backlit

## Automatic Temperature Compensation

0°C~50°C

## Environment

Temperature:0°C~50°C; Humidity:≤95%

## Power Supply

4\*1.5V(AG13)

## Dimensions

170mm\*33mm\*20mm

## Weight

69g