

Digital Turbidity Meter



Product Categories: [Analytical Lab Instruments](#), [Calorimeters](#)

Product Tags: [chemistry lab instruments](#)

Product Page: <https://www.labappara.com/product/digital-turbidity-meter-2/>

Product Description

DIGITAL TURBIDITY METER - Chemistry Equipment

Three & half **LED** display.

Range up to 1000 NTU/JTU. Resolution 1 NTU/JTU.

Readout

3 ½ digit LED display

Range

0 - 200 NTU

0 - 2000 NTU

Resolution

0.1 NTU in 200 range

1 NTU in 2000 range

Accuracy

± 2% FSR

Power required

230 V AC ± 10% 50 Hz with inbuilt stabilizer

Warming time

5 minutes

Weight

1 Kg (approx.)

Body

ABS

Dimension

95 mm (H) x 225 mm (W) x 215 mm (L)

Accessories

Beaker set of 2 pieces, 2 bottles of standard solution (No. I and II) and Dust proof cover