

Digital Water And Soil Analysis Kit



Product Categories: [Analytical Lab Instruments](#)

Product Tags: [Analytical Laboratory Equipment](#)

Product Page:

<https://www.labappara.com/product/digital-water-soil-analysis-kit-2/>

Product Description

Digital Water And Soil Analysis Kit - Analytical Laboratory Equipment

Three and half digit LCD display covering Six parameter pH, Conductivity, TDS, ORP, Salinity and Temperature.

Battery cum mains operated, complete with accessories, fitted in briefcase.

Range : 0-14 pH

Resolution : 0.01pH

Accuracy : + 0.01 pH (relative)

Temp. Comp. : Automatic with temp. Probe

Buffer : Manuals with keys

Recognition : 4.00 , 7.00 , 9.20 (Automatic)

Slope Adjustment : 80% to 120%, by default 100%

Calibration : Auto/Manual 2-Point, 3-Point