

## Diode Valve Characteristics Apparatus



**Product Categories:** [E/M Apparatus & Valve Experiments](#), [E/M Apparatus & Valve Experiments](#), [Electronics](#), [Engineering Equipment](#)

**Product Page:**

<https://www.labappara.com/product/diode-valve-characteristics-apparatus/>

### Product Description

#### Diode Valve Characteristics Apparatus

**Features:**

Instrument comprises of DC Regulated Power Supply 0-250VDC/30mA, Variable AC Power Supply 6.3 VAC/1Amp, two round meters for voltage & current measurement.