

Dispensing Balance



Product Categories: [Balance & Scale](#), [Chemistry Equipments](#)

Product Tags: [chemistry lab instruments](#)

Product Page: <https://www.labappara.com/product/dispensing-balance/>

Product Description

Dispensing Balance - Laboratory Equipment

Dispensing Balance is a 2 pan torsion type balance which utilizes both internal and external weights. It requires the use of external weights for measurements

Capacity:- 50 gm, 100 gm Feature:-

Double pan packed in wooden box without weight

[Chemicals](#) or small objects up to 100 gm can be weighed to within 15 mg.

The brass beam support is mounted on the wooden base.

Equipped with a drawer for holding weights and other equipment's.