

Display Interface



Product Categories: [Electronics](#), [Engineering Equipment](#), [Interfacing Cards](#), [Interfacing Cards](#), [Microprocessor & Micro-controller Lab](#), [Microprocessor & Micro-controller Lab](#)

Product Page: <https://www.labappara.com/product/display-interface/>

Product Description

Display Interface Specification: 6 No of Seven Segment

Matrix Key Interface: 18 No of snap on keys

Simultaneous Keyboard and Display operation.

3 input mode such as scanned keyboard mode, scanned sensor mode and strobed input entry mode

6 character Display

8 character keyboard FIFO

Right or Left entry 6 Byte Display RAM

Programmable Scan Timing

2 Key lockout or N-key Roll-over with contact debounce.

Auto Increment facility for easy programming

50 pin FRC connector to interface with 8085/8086/8051 Trainer Kits

Features:

Accurate dimension

Longer service life

Light weight

Maintenance free