

Dissecting Microscope(DM-01)



Product Categories: [Microscope](#), [Microscopes](#)

Product Tags: [biology lab equipment](#)

Product Page: <https://www.labappara.com/product/dissecting-microscope-2/>

Product Description

Dissecting Microscope(DM-01) - Scientific Equipment

Dissecting Microscope has a heavy round base with precisely designed body. Sensitive focusing is done by rack and pinion arrangement. Revolving arm provided for moving magnifying lens over full stage area. It is easily replaceable by a special bull's eye lens (83 mm dia.) to increase the field of vision. A piano-concave mirror reflects sufficient day light and artificial light on the object. Detachable hand rests on both sides of stage provided for a study support during dissection, provided with 10 x & 20 x eye pieces and special bull's eye lens (83 mm dia.), packed in wooden carrying cabinet.