

Dissecting Microscope



Product Categories: [Microscope](#), [Microscopes](#)

Product Tags: [biology lab equipment](#)

Product Page: <https://www.labappara.com/product/dissecting-microscope-4-2/>

Product Description

Dissecting Microscope(DM 04) - Scientific Equipment

Dissecting Microscope with bull lens. Housed in thick plywood box.

Widely used for study of small organisms and for Dissection work.

Heavy round C.I. base for greater stability.

All brass chrome plated parts like pipe, knobs & lens holding arm.

Wide glass stage having two stage clips & detachable Hand Rests on both sides.

Focusing by Rack & Pinion mechanism.

Eye Pieces 10 x & 20 x made out of OPTICAL GLASS are mounted on joint revolving arm

allowing the whole stage area to be scanned. Magnification 10 x and 20 x.

Illumination by Plano-Concave reflector.

Packed in a sturdy & beautiful Sun-mica clad plywood box with dust cover & cleaning cloth.

Optical Combination:

Eye Piece : 10 X & 20 X Eye Piece

Object stage size : 110 mm x 115 mm

Features:

Superior functionality

Long service life

Easy maintenance

User friendly

Crystal clear lenses

Light weight