

## Dissecting Microscope(DM- 03)



**Product Categories:** [Microscope](#), [Microscopes](#)

**Product Tags:** [biology lab equipment](#)

**Product Page:** <https://www.labappara.com/product/dissecting-microscope-4/>

### Product Description

Dissecting Microscope(DM- 03) - Scientific Equipment

Dissecting Microscope with round base, brass knobs, brass eyepiece with aluminium fitting eyepiece 10 x, 20 x . Well packed in Thermocol box.

#### **Optical Combination:**

Eye Piece : 10 X & 20 X Eye Piece

Object stage size : 110 mm x 115 mm

#### **Features:**

Superior functionality

Long service life

Easy maintenance

User friendly

Crystal clear lenses

Light weight