

Distillation Anti Splash Adapter



Product Categories: [Distilling Apparatus Borosilicate Glass](#), [Laboratory Glassware](#)

Product Tags: [glassware](#), [Laboratory Equipment](#)

Product Page: <https://www.labappara.com/product/distillation-anti-splash-adapter/>

Product Description

Distillation Anti Splash Adapter - Laboratory Glassware Equipment

Designed with inlet and outlet neck, the offered adapter is extensively used for connecting two or more tubes or apparatus together. It is also known as interconnecting component. Accessible in different sizes, this Splash Head Adapter Straight & Bend can be availed from us at feasible prices.

A splash head for use in distillation assemblies in which the condenser is assembled in the vertical mode.

Lower cone accepts flask, upper cone accepts condenser.

Splash head reduces entertainment of raw liquid into the vapor stream.

Features:

Excellent strength

Leakage proof

Transparency

Fine finish