

Distillation Columns



Product Categories: [Distilling Apparatus Borosilicate Glass](#), [Laboratory Glassware](#)

Product Tags: [glassware](#), [Laboratory Equipment](#)

Product Page: <https://www.labappara.com/product/distillation-columns/>

Product Description

Distillation Columns - Laboratory Glassware Equipment

Columns acc. To hempel for packing like raschig rings or glass beads, with ns socket and ns cone, without glass jacket, material: borosilicate glass 3. 3

Distillation is a process that can be used to separate a pure liquid from a mixture of liquids.

It works when the liquids have different boiling points.

Distillation is commonly used to separate ethanol - the alcohol in alcoholic drinks - from water. The mixture is heated in a flask.