

## Distillation Head



**Product Categories:** [Distilling Apparatus Borosilicate Glass](#), [Laboratory Glassware](#)

**Product Tags:** [glassware](#)

**Product Page:** <https://www.labappara.com/product/distillation-head/>

### Product Description

#### Distillation Head - Laboratory Glassware Equipment

Distillation is the separation of a mixture into its component parts, or fractions, separating chemical compounds by their boiling point by heating them to a temperature at which one or more fractions of the compound will vaporize. This head is mainly used to quickly fit or join different apparatus together. Tested by quality controllers to ensure its durability. Distilling Head made of Borosilicate 3.3 glass which is resistant to heat & almost many chemicals, suitable to be used in Distillation Assemblies for connecting Flasks to condensers with side cone at 75° angled to the center.

#### Features:

Durable

Light weight

Easy to connect

## Flawless finish