

# Distillation Vertical



**Product Categories:** [Distilling Apparatus Borosilicate Glass](#), [Laboratory Glassware](#)

**Product Tags:** [glassware](#)

**Product Page:** <https://www.labappara.com/product/distillation-vertical/>

## Product Description

Distillation Vertical - Laboratory Glassware Equipment

Distillation is a process that can be used to separate a pure liquid from a mixture of liquids.

It works when the liquids have different boiling points.

Distillation is commonly used to separate ethanol - the alcohol in alcoholic drinks - from water.

The mixture is heated in a flask.