

Dressing Trolley



Product Categories: [Hospital Equipment & Furniture](#)

Product Tags: [hospital](#), [medical](#)

Product Page: <https://www.labappara.com/product/dressing-trolley/>

Product Description

Dressing Trolley - Hospital Equipment

This trolley finds its applications in nursing homes and hospitals for keeping dressing materials and surgical tools. Under expert's observation, the offered trolley is manufactured using top-of-the-line stainless steel.

Frame Work Constructed With Stainless Steel Material.

Top Of The Trolley Covered With The Stainless Steel Rod From Three Sides And The Bottom Side Was Covered With SS Rod With Four Sides.

Two Stainless Steel Shelves at Bottom and Top.

Overall Size - 75 L X 45 W X 90 H cm