

Dry & Wet Thermometer (Duplex Box)



Product Categories: [Chemistry Equipments](#), [Thermometer & Hygrometer](#)

Product Tags: [chemistry lab instruments](#), [physics lab instruments](#)

Product Page: <https://www.labappara.com/product/dry-wet-thermometer-duplex-box/>

Product Description

Dry & Wet Thermometer (Duplex Box) - Scientific Equipment

We are offering Wet thermometer, Dry Thermometer.

We are engaged in manufacturing and exporting a wide range of wet and dry thermometers.

These are provided with mercury filled imported glass and white back thermometers, fixed on metal scale.

Also, these thermometers are equipped with superior quality reinforced polypropylene plastic case, and a lid & bulb protector.

These are also provided with removable spare wicks and mounting screws & plastic water cistern and boxed with hygrometer tables.

Features:

Superb functionality

Effective

Reliable
Easy to use
Long service life