

E/M By Helical's Method



Product Categories: [E/M Apparatus & Valve Experiments](#), [E/M Apparatus & Valve Experiments](#), [Electronics](#), [Engineering Equipment](#)

Product Page: <https://www.labappara.com/product/e-m-by-helicals-method/>

Product Description

E/M By Helical's Method

Features :

Kit Comprises of High Voltage Power Supply Unit with intensity, focus X, Y, deflection & Solenoid current controls. Two meters are provided in Power Supply for Acceleration Voltage control & for solenoid current controls. One Solenoid Unit Comprises of one Solenoid mounted on Wooden stand & CRT is kept in the solenoid pipe.