

## E.S.R Tube



**Product Categories:** [Pathology Lab](#)

**Product Tags:** [Physics equipment](#)

**Product Page:** <https://www.labappara.com/product/e-s-r-tube/>

### Product Description

#### Erythrocyte Sedimentation Rate (E.S.R) Tube-Physics Equipment

The Erythrocyte Sedimentation Rate (ESR) is a common hematological test for nonspecific detection of inflammation that may be caused by infection, some cancers and certain autoimmune diseases.

These tubes are to be used with various automated Erythrocyte sedimentation rate analyzers.

After adding the whole blood in to the tube mix thoroughly for 7 - 8 times.

It can be defined as the rate at which Red Blood Cells sediment in a period of one hour.