

Electric Bell



Product Categories: [Educational equipments](#), [Electromagnetism](#), [Physics Equipments](#)

Product Tags: [physics lab instruments](#)

Product Page: <https://www.labappara.com/product/electric-bell/>

Product Description

Electric Bell - Physics Working Model

For Demonstration all parts easily visible, mounted on sun mica top base workable 4-6 volt DC.

Specifications:

Power Source: Electric

Voltage: 4 - 6 Volt

Material Used: SS, Copper, Wood

Power Source

Electric

Voltage

6V

Material Used

SS, Copper, Wood

Application

Physics Lab