

Electrical ICU Bed



Product Categories: [Hospital Equipment & Furniture](#)

Product Tags: [hospital](#), [medical](#)

Product Page: <https://www.labappara.com/product/electrical-icu-bed/>

Product Description

Electrical ICU Bed - Hospital Equipment & Furniture

Framework made of rectangular CRCA tubes section top made of perforated M.S. sheet.

Backrest, knee rest, and hi-lo positions maneuvered by separate screw from foot end.

Collapsible / swing away safety side railing

Telescopic IV rod with four locations

Bed mounted on 125 mm diameter castor two with brakes.

Rubber buffers at four corners.

Pre-treated and epoxy power coated finish.

Overall Size: 218 (L) x 98 (W) x 55-80 (H) cm