

Electromagnet U- Shaped



Product Categories: [Educational equipments](#), [Electromagnetism](#), [Physics Equipments](#)

Product Tags: [physics lab instruments](#)

Product Page: <https://www.labappara.com/product/electromagnet-u-shaped/>

Product Description

Electromagnet U- Shaped - Physics Equipment

Electromagnet U shape has two coils of copper and inserted in U shape iron rod. Comprising iron U-shaped core of circular section with armature and carrying hook. Armatures on both the arms wound with enameled copper wire and are joined in series. Operates on 3-6 Volts D.C. Overall length 100 mm.

Available in size - 25 mm, 50 mm, 75 mm, 100 mm, 125 mm, 150 mm