

Electronic Balance Analytical



Product Categories: [Balance & Scale](#), [Chemistry Equipments](#)

Product Tags: [biology lab equipment](#), [chemistry lab instruments](#), [physics lab instruments](#), [science lab equipment](#)

Product Page: <https://www.labappara.com/product/electronic-balance-analytical/>

Product Description

Electronic Balance Analytical - Scientific Equipment]

Electronic analytical scales measure the force needed to counter the mass being measured rather than using actual masses.

Feature - Highly accurate - Calibrated - Widely used in laboratories - Consistent performance - Sturdiness - Competitive rate - Low maintenance