

Expansion Adapters



Product Categories: [Adapters Borosilicate Glass](#), [Laboratory Glassware](#)

Product Tags: [glassware](#), [Laboratory Equipments](#)

Product Page: <https://www.labappara.com/product/expansion-adapters/>

Product Description

Expansion Adapters - Laboratory Glassware Equipment

We hold all the aces in developing the best Reduction Expansion Adapters. Manufactured by making use of superior raw-material and experienced work force, these products are widely appreciated for their durability and high efficiency. In order to meet diverse demands of our valued clients, we offer these adapters in different sizes that can be customized as per their requirements.

All expansion adapters are made from borosilicate glass. They are useful for the connection of dissimilar sized ground glass joints. Internchangeable large socket at the top and interchangeable small cone at the base.

Features :

Long lasting

Sturdiness

Effectiveness

Easy to use

Fine-finish

Rigid

Smooth threads

Elegantly Designed

No wear and tear

Super strength

Super strengthen

Great efficiency

Handles high pressure

Long lasting property

Available in socket size - 19/26, 24/29, 29/32, 34/35, 40/38