

Extraction Apparatus Borosilicate Glass



Product Categories: [Extraction Apparatus Borosilicate Glass](#), [Laboratory Equipments](#), [Laboratory Glassware](#)

Product Tags: [glassware](#), [Laboratory Equipments](#)

Product Page:

<https://www.labappara.com/product/extraction-apparatus-borosilicate-glass/>

Product Description

Extraction Apparatus Borosilicate Glass- Laboratory Glassware Equipment Laboratory", consists of Flask, Extractor and Condenser.

Extraction Apparatus is commonly used to separate the compounds from insoluble substances. extraction glass is the most suitable laboratory apparatus. The quality controlled calibration offers reliability in results due to complete accuracy.

Available Extraction Capacity - 20 ml, 40 ml, 60 ml, 100 ml, 100 ml, 200 ml, 400 ml, 600 ml, 1000 ml, 2000 ml, 3000 ml.

Extractor

Capacity

ml

Extractor

Socket

Extractor

Cone

Flask

Capacity

R.B./F.B.

20

24/29

19/26

50

40

29/32

24/29

100

60

34/35

24/29

150

100

34/35

24/29

250

100

40/38

24/29

250

200
50/42
24/29
500

400
50/42
24/29
1000

600
50/42
24/29
2000

1000
55/44
34/35
3000

2000
55/44
34/35
5000

3000
55/44
34/35
10000