

FET Characteristics Apparatus



Product Categories: [Characteristics And Application](#), [Electronics](#)

Product Page: <https://www.labappara.com/product/fet-characteristics-apparatus/>

Product Description

FET Characteristics Apparatus

Features:

Two moving coil voltmeters MO 65 mounted on the front panel to measure voltage. The connections for each meter is brought out at 4mm sockets

One moving coil ammeters MO 65 mounted on the front panel to measure Drain Current. The connections for each meter is brought out at 4mm sockets

Two DC Regulated Continuously Variable Power Supplies 0-15 Volts DC / 150 mA and 0-5 Volts DC / 150 mA

One FET (Field Effect Transistor) mounted behind the panel and the connection for Drain Source and Gate are brought out at 4mm sockets.

One ON/OFF Switch with jewel light is provided on the front panel.