

## Filter Apparatus



**Product Categories:** [Laboratory Glassware](#), [Sintered Ware Borosilicate Glass](#)

**Product Tags:** [glassware](#)

**Product Page:** <https://www.labappara.com/product/filter-apparatus/>

### Product Description

Filter Apparatus - Laboratory Glassware Equipment

Filter Apparatus Complete, Filter Flask with Socket and buchner

Borosilicate glass with fritted glass filter base.

Base design has integral vacuum connection located above filter drip to prevent contamination of vacuum line with droplets.

funnel with Sintered Disc of Porosity G-1 or G-2 or G-3 or G-4, having cone

#### **Features:**

Glass filtration system for use with 47 mm diameter membranes

Borosilicate glass sintered disc to support membrane

Funnel 300 ml capacity

Support base with sintered disc and vacuum connection

Fully autoclave

Cat. No.	Flask Cap.	Buchner Funnel
ml	Cap. ml	
327/1	250	80
327/2	500	200
327/3	1000	500
327/4	2000	1000