

# Filter Funnel Buchner



**Product Categories:** [Funnels Borosilicate Glass](#), [Laboratory Glassware](#)

**Product Tags:** [glassware](#), [Laboratory Equipment](#)

**Product Page:** <https://www.labappara.com/product/filter-funnel-buchner/>

## Product Description

Filter Funnel Buchner - Laboratory Glassware Equipment

Vacuum filtration is a technique for separating a solid product from a solvent or liquid reaction mixture.

The mixture of solid and liquid is poured through a filter paper in a Buchner funnel.

The solid is trapped by the filter and the liquid is drawn through the funnel into the flask below, by a vacuum.

**Features:**

Dimensional accuracy

Equipped with glass sintered disc

High resistance to chemical attack