


Filter Funnel Buchner



Product Categories: [Funnels Borosilicate Glass](#), [Laboratory Glassware](#)

Product Tags: [glassware](#), [Laboratory Equipment](#)

Product Page: <https://www.labappara.com/product/filter-funnel-buchner/>

Product Description

Filter Funnel Buchner - Laboratory Glassware Equipment

Vacuum filtration is a technique for separating a solid product from a solvent or liquid reaction mixture.

The mixture of solid and liquid is poured through a filter paper in a Buchner funnel.

The solid is trapped by the filter and the liquid is drawn through the funnel into the flask below, by a vacuum.

Features:

Dimensional accuracy

Equipped with glass sintered disc

High resistance to chemical attack