


Filter Funnel



Product Categories: [Funnels Borosilicate Glass](#), [Laboratory Glassware](#)

Product Tags: [glassware](#), [Laboratory Equipment](#)

Product Page: <https://www.labappara.com/product/filter-funnel/>

Product Description

Filter Funnel - Laboratory Glassware Equipment

A filter funnel is a laboratory funnel used for separating solids from liquids via the laboratory process of filtering. In order to achieve this, a cone-like shaped piece of filter paper is usually folded into a cone and placed within the funnel.

The suspension of solid and liquid is then poured through the funnel.

Features:

High strength

Precise Design

Easy to handle

Dia (cms) 5, 6.25, 7.5, 10.15

Accurate measurement

Transparent