

Fish Tail Burner



Product Categories: [Chemistry Equipments](#), [Lab Supplies](#)

Product Tags: [chemistry lab instruments](#), [science lab equipment](#)

Product Page: <https://www.labappara.com/product/fish-tail-burner/>

Product Description

Fish Tail Burner - Scientific Instrument

A Fish tail burner in which two jets of gas issuing from two holes inclined toward each other impinge and form a fan-shaped flame. Chrome plated brass tube with non-adjustable jet. On base with ruffled gas inlet tube. One spare gas nipple is provided for use with LPG gas.

Feature

On cast iron base.

Wide & stable base provide

Packed in box