

Fisher Clamp



Product Categories: [Clamps, Holders Racks Stands](#), [Laboratory Plastic Ware](#)

Product Tags: [biology lab equipment](#), [chemistry lab instruments](#), [Manufacturer of Laboratory equipment](#), [math lab instruments](#), [physics lab instruments](#), [science lab equipment](#)

Product Page: <https://www.labappara.com/product/fisher-clamp/>

Product Description

Fisher Clamp - Scientific Laboratory Equipment

The Fisher Clamps are molded in Polypropylene making it completely corrosion-proof. The spring-loaded jack in the middle provides a vice like grip to the Burette and prevents it from slipping.

Material : Polypropylene

Features:

Excellent strength

Flawless finish

Low maintenance

Particulars	Packing
Fisher Clamp - Single	12 Pcs
Fisher Clamp - Double	6 Pcs