

Fixed Frequency Sine Wave Oscillator (1 KHz)



Product Categories: [Electronics](#), [Engineering Equipment](#), [Function Generators](#), [Function Generators](#)
Product Page: <https://www.labappara.com/product/fixed-frequency-sine-wave-oscillator-1-khz/>

Product Description

Fixed Frequency Sine Wave Oscillator (1 KHz)

Using operational amplifier, Wein Bridge oscillator for better stability.

It provides pure sine wave 0-10 volt R.M.S. continuously variable output, giving fixed frequency 1KHz.

This bridge oscillator is used in almost all the four arm bridges.