

FLAME PHOTOMETER DOUBLE DISPLAY



Product Categories: [Analytical Lab Instruments](#), [Spectrophotometer](#)

Product Tags: [Analytical Lab Equipment](#), [Spectrophotometer](#)

Product Page:

<https://www.labappara.com/product/flame-photometer-double-display/>

Product Description

Flame Photometer- Double Display- Analytical Lab Equipment

2 Filter Na, k(LI& Ca Optional)

Range from 0 to 100ppm

Narrow band interference glass filters

LPG and dry oil free air

Double display 2 ½ digit bright red led Display.

Measure one element at a time.

Available in Single Display