

# Flashing & Quenching Of Neon Bulb



**Product Categories:** [Characteristics And Application](#), [Characteristics And Application](#), [Electronics](#), [Engineering Equipment](#)

**Product Page:** <https://www.labappara.com/product/flashing-quenching-of-neon-bulb/>

## Product Description

Flashing & Quenching Of Neon Bulb

**Features :**

Instrument comprises of High Voltage DC Power Supply 0-250 VDC/30mA, 3 Capacitors (Paper Type) connected behind the front panel, can be selected using toggle switches & one Neon Lamp.