

Flask Stand



Product Categories: [Clamps, Holders Racks Stands](#), [Laboratory Plastic Ware](#)

Product Tags: [biology lab equipment](#), [chemistry lab instruments](#)

Product Page: <https://www.labappara.com/product/flask-stand/>

Product Description

Flask Stand - Scientific Laboratory Plastic ware

Molded in Polypropylene, the Flask Stand gives a stair-like top view. Every step provides excellent support to the round bottom flask at every point of Contact. This Stackable Flask Stand can also be autoclaved.

Particulars	Packing
-------------	---------

Flask Stand	6 Pcs
-------------	-------