


Four Finger Clamp



Product Categories: [Chemistry Equipments](#), [Lab Supplies](#)

Product Tags: [chemistry lab instruments](#), [physics lab instruments](#), [science lab equipment](#)

Product Page: <https://www.labappara.com/product/four-finger-clamp/>

Product Description

Four Finger Clamp - Scientific Instrument

Four finger clamps are used to hold larger pieces of glassware and odd-shaped objects.

Feature -

Brass oxidised with boss head.

Provide with 7" mild steel rod C.P.

Clamp is staple painted, jaws are cork lined.

Made of Aluminum Alloy

Made of brass casted and C.P.

Jaws are Made aluminium casted material.

Jaws are crossed and duly covered with thick rubber sheet.

Able to catch up to 75/125 mm dia objects

150/200/250mm long chrome plated rod

Duly powder coated finishing.

Well designed

Best grade components

Latest technology