

Gas Generator Borosilicate Glass



Product Categories: [Gas Generator Borosilicate Glass](#), [Laboratory Glassware](#)

Product Tags: [glassware](#), [Laboratory Equipment](#)

Product Page:

<https://www.labappara.com/product/gas-generator-borosilicate-glass/>

Product Description

Gas Generator Borosilicate Glass - Laboratory Glassware Equipment

It is made up of borosilicate glass with ultra-modern technique to ensure its optimum quality and reliability. This is an effective and rapid device designed to generate controllable amounts of common gases like carbon dioxide, hydrogen sulphide, etc. This apparatus can be swiftly changed & assembled and is easy to operate. Equipped with thistle funnel and short delivery tube, the Gas Generator. We are different size and shapes and trading Gas Generator Kipp's Borosilicate Glass. It is used for executing many chemical reactions for determination of specific metal ions.

Features:

Resistant to heat

Fine finish

High strength

Capacity: 500 ml