

Gas Taps



Product Categories: [Chemistry Equipments](#), [Lab Supplies](#)

Product Tags: [chemistry lab instruments](#), [science lab equipment](#)

Product Page: <https://www.labappara.com/product/gas-taps/>

Product Description

Gas Taps- Scientific Equipment

We are engaged in offering Gas Tap, which are made of brass and are used in hospitals, colleges, laboratories and research institutes. All brass black powder coated, Heavy Leak proof stop cock. with adopter.

Feature

Brass C.P guaranteed leak proof in card box.

Elegant design

Scratch resistant

Chemical paint

Chrome plated

Standard 3 mm

Inlet 1/4 inches

54mm diamete

Provide with one way, two way, three way, four way.