

Glass Cuvette with PTFE Cover



Product Categories: [Laboratory Glassware](#), [Quartz Ware](#)

Product Tags: [glassware](#)

Product Page: <https://www.labappara.com/product/glass-cuvette-ptfe-cover/>

Product Description

Glass Cuvette with PTFE Cover - Laboratory Glassware Equipment

All scientists know that their laboratory must be equipped with a few standard 10 mm cuvettes. Precision Cells Type 1 standard 10 mm cell has become our best seller for spectroscopy and with many different types of materials available it's easy to find the exact cuvette you need.

Features:

Highly accurate results, only need to test sample one time

Accurate repeat-ability, Precision Cells cuvettes are made with only the highest quality quartz and glass materials

Hand made, ensuring the highest quality optical polish with no defects.

Every cell is thoroughly inspected before leaving the factory.

Precision Cells cuvettes are guaranteed to give you the best optical transmission.

Fire fused to ensure the best bonding.

Three premium materials available: (Optical Glass, Pyrex Glass (10mm only), UV Quartz, IR Quartz)

Lightpaths from 1 mm up to 100 mm available.

Two windows optically polished.

Includes PTFE cover.

Place your order on-line using the chart below.

Recommended accessories:

Spacer for adapting 1, 2 and 5 mm versions to a standard 10 mm cuvette holder:

A25 Cuvette Spacer

Spring-based mount for fitting any versions less than 5mm to a standard 10

mm cuvette holder: A2o Series: Cuvette Mount

Cuvette Holder for cells with 1-100 mm path-length: CHR1/x

Mixer for use with 10 mm cells: P69 Mixer

Magnetic stirrer for use with 10 mm cells: P73 Stirrer

Cleaning swabs: SwabMFC

Spare PTFE covers

10 mm path, 3.5 ml in pair.