

Glass Prism



Product Categories: [Educational equipments](#), [Glass Prism](#), [Light & Optics](#), [Optical Instruments](#), [Physics Equipments](#)

Product Tags: [physics lab instruments](#)

Product Page: <https://www.labappara.com/product/glass-prism/>

Product Description

Glass Prism - Physics Optical Equipment

Glass prism reflection Prism's equilaterals have and also used to disperse.

This kind of Glass Prism can reflect seven-color sunlight.

Great prism for teaching can be also a toy for kids made of high quality optical glass.

Size (a) 38x 38mm

(b) 50x 50mm