

## H.B Meter



**Product Categories:** [Pathology Lab](#)

**Product Tags:** [Laboratory Equipment](#)

**Product Page:** <https://www.labappara.com/product/h-b-meter/>

### Product Description

H.B Meter - Laboratory Equipment

Comparator Holder: -1 Nos.

It is black in colour. The back of the holder is fitted with milky colour screen made of acrylic. In front of the holder there are three windows. The two side windows are fitted with colour comparison. The size of all the three windows are Height 3 cms. (Minimum) Width 0.5 cms. (Minimum)

Square H.B. Tube: -2 Nos.

Square H.B. Tube is graduated on both the side for the measurement of haemoglobin in respect of percentage and gram.

H.B. Pipette -1 Nos.

20 microlite H.B. Pipette with latex tube and mouth piece

Amber coloured glass bottle -1 Nos.

Glass Stirrer -1 Nos.

Dropper with teat -1 Nos.

Cleaning Brush -1 Nos.

Spare latex tubing one meter -1 Nos.

Hemoglobin is the protein molecule in red blood cells that carries oxygen from the lungs to the body's tissues and returns carbon dioxide from the tissues back to the lungs.

Hemoglobin is made up of four protein molecules ( globulin chains ) that are connected together.

**Feature -**

Hemoglobin is consists of Regular quality.

Hemoglobin is consists of superior quality in German made.