


# Hamming Code Trainer



**Product Categories:** [Advance Communication Lab](#), [Engineering Equipment](#)  
**Product Page:** <https://www.labappara.com/product/hamming-code-trainer/>

## Product Description

Hamming Code Trainer

### Features:

- Rugged design
- Light weight
- Enhanced performance

### Other Details:

- Onboard transmitter clock and Sync (Fast and Slow mode)
- On Board Modulating Digital Data signal generator to generate any binary input word with Bit clock & Word Clock. Word Length: 7 Bits
- On Board check bits: Three bits and Four bits data
- On Board Noise inserting circuit through four bit error generator
- Onboard seven bit shift register
- To generate digital data word with four bit data & three parity check bits