

Haemocytometer



Product Categories: [Pathology Lab](#)

Product Tags: [physics Equipments](#)

Product Page: <https://www.labappara.com/product/haemocytometer/>

Product Description

Haemocytometer

The (Haemocytometer or counting chamber) is a specimen slide which is used to determine the concentration of cells in a liquid sample. It consists of 2 pipettes, one for red blood corpuscles and the other for white blood corpuscles, and a Counting Chamber with double Neubauer ruling, in a case. The Counting Chamber is made from a single piece of glass with a ruling engraved directly on the glass and with wide protecting and transverse moats. Above and below the ruled area is a matte surface, which provides a mean for easy focusing. A location line facilitates immediate finding of the rulings in the microscope field without disturbing the preparation or damaging the objective by touching the cover glass.

Feature -

Complete with German chamber & W.B.C & R.B.C Pipette original.

With china chamber & W.B.C & R.B.C Pipette regular quality