

Hirsch Funnel



Product Categories: [Educational equipments](#), [Laboratory Porcelain Ware](#)

Product Tags: [biology lab equipment](#), [chemistry lab instruments](#), [physics lab instruments](#)

Product Page: <https://www.labappara.com/product/hirsch-funnel/>

Product Description

Hirsch Funnel - Scientific Laboratory Equipment

Superior quality chemical- porcelain Hirsch Funnels are designed with conical sides, a fixed perforated plate, and a tapered stem. Glazed inside and outside, aside from the rim, they are highly chemical resistant and autoclave. Generally used to filter particulates from solutions.

Diameter :- 25 mm, 35 mm, 50 mm, 60 mm, 75 mm, 100 mm.

To prolong the life of porcelain products:

Follow gradual heating/cooling rates

Use an oven or hot plate and an intermediate step when quicker heating/cooling rates are required

Avoid contact of heated ware with a cold surface

Carefully inspect lab ware prior to each use and do not use any product that appears defective