

Hoffman Screw Clip



Product Categories: [Chemistry Equipments](#), [General Chemistry Labware](#)

Product Tags: [chemistry lab instruments](#), [physics lab instruments](#), [science lab equipment](#)

Product Page: <https://www.labappara.com/product/hoffman-screw-clip/>

Product Description

Hoffman Screw Clip(iron) - Scientific Equipment

Specifications:

Material: Copper, Iron

Screwed, hinged side bar, for use with tubes upto 20 mm outer diameter

Plated metal

Plated brass

Features:

Excellent corrosion resistance

High temperature withstanding capabilities

Smooth pore and surface

Maintenance free

Increased durability

Catalytic, inertness

Physiological inertness

High tensile strength

Chemical resistant