

Horizontal Pipette Stand



Product Categories: [Laboratory Plastic Ware](#), [Pippetes Stand](#)

Product Tags: [Laboratory Equipments](#), [Laboratory plastic ware](#)

Product Page: <https://www.labappara.com/product/horizontal-pipette-stand/>

Product Description

Horizontal Pipette Stand - Scientific Laboratory Equipment

Pipette Stand can hold 12 pipettes, 6 on each side and is detachable and autoclavable. Moulded in Polypropylene. Horizontal Pipette Stand has two sides that are attached with each other with three metallic rods coated with Polypropylene. For storage purpose stand can be easily dismantled.

Features:

High strength

Perfect finish

Elevated durability

Particulars

Packing

For 12 Pipettes

6 Pcs