

## Horse Shoe Magnet (Alnico)



**Product Categories:** [Educational equipments](#), [Electromagnetism](#), [Physics Equipments](#)

**Product Tags:** [physics lab instruments](#)

**Product Page:**

<https://www.labappara.com/product/horse-shoe-magnet-alnico-type/>

### Product Description

Horseshoe Magnet is Highly effective when high magnetic strength is required. These magnets will hold their force for a long time. Distance between the poles is 20 mm.

#### **Features:**

Robust Construction

Perfect Finish

Enhanced functional life

Fine quality

Easy to maintain

## Durability

Available in size - 25 mm, 50 mm, 75 mm, 100 mm, 125 mm, 150 mm