

Horse Shoe Magnet



Product Categories: [Educational equipments](#), [Electromagnetism](#), [Physics Equipments](#)

Product Tags: [physics lab instruments](#)

Product Page: <https://www.labappara.com/product/horse-shoe-magnet/>

Product Description

Horse Shoe Magnet - Physics Equipment

A **horseshoe magnet** is a **magnet** made in the shape of **horseshoe**. The **magnet** has two **magnetic** poles close together. This shape creates a strong **magnetic** field between the poles. Our firm is providing a best class **Horseshoe Magnet** to the prestigious customers. A **U-shaped magnet** derives its name from its **shape** and has both a north and a south pole located in the same plane at the open end of the **magnet**. This type of **magnet** is also called a **horseshoe magnet**. These are of steel with keepers, highly polished and magnetized)

Features:

- Easy to use - High efficiency - Lightweight - Trouble-free performance.

Available in size :- 25mm,50mm ,75mm , 100mm ,125mm, 150mm