

Hospital Examination Couch Technical Specification



Product Categories: [Hospital Equipment & Furniture](#)

Product Tags: [hospital](#), [medical](#)

Product Page:

<https://www.labappara.com/product/hospital-examination-couch-technical-specification/>

Product Description

Hospital Examination Couch Technical Specification - Hospital Equipment

Constructed from CRCA steel sheet.

Manually Adjustable Head Section.

Two individual cabinets with attached drawers

Epoxy power Coated.

Overall Size - 1830 L X 60 W X 810 H